





TGFBR1 Antibody

Product Code	CSB-PA004276
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36897
Immunogen	Synthesized peptide derived from Human TGFβ RI around the non-phosphorylation site of S165.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TGFBR1; ALK5; SKR4; TGF-beta receptor type-1; TGFR-1; Activin A receptor type II-like protein kinase of 53kD; Activin receptor-like kinase 5; ALK-5; ALK5; Serine/threonine-protein kinase receptor R4; SKR4; TGF-beta type I receptor; Transfor
Immunogen Species	Homo sapiens (Human)
Target Names	TGFBR1

Image



Western Blot analysis of HepG2 cells using TGFβ RI Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using
TGFβ RI Polyclonal Antibody