





TGFBR1 Antibody

Product Code	CSB-PA004276
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36897
Immunogen	Synthesized peptide derived from Human TGFβ RI around the non-phosphorylation site of S165.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGFBR1

Image



Western Blot analysis of HepG2 cells using TGFβ RI Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using TGFβ RI Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.