





TGFBR2 Antibody

Product Code	CSB-PA004277
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37173
Immunogen	Synthesized peptide derived from Human TGFβ RII around the non-phosphorylation site of S225.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TGFBR2; TGF-beta receptor type-2; TGFR-2; TGF-beta type II receptor; Transforming growth factor-beta receptor type II; TGF-beta receptor type II; TbetaR-II
Immunogen Species	Homo sapiens (Human)
Target Names	TGFBR2

Image



Western Blot analysis of NIH-3T3 cells using TGFβ RII Polyclonal Antibody



Western blot analysis of 3T3 lysis using TGFβ RII antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.