




TH Antibody

Product Code	CSB-PA004286
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07101
Immunogen	Synthesized peptide derived from Mouse TH around the non-phosphorylation site of S8.
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TH; TYH; Tyrosine 3-monooxygenase; Tyrosine 3-hydroxylase; TH
Immunogen Species	Rattus norvegicus (Rat)
Target Names	TH

Image



Western Blot analysis of COLO205 cells using TH Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.