




TH Antibody

Product Code	CSB-PA004286
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07101
Immunogen	Synthesized peptide derived from Mouse TH around the non-phosphorylation site of S8.
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TH; TYH; Tyrosine 3-monooxygenase; Tyrosine 3-hydroxylase; TH
Immunogen Species	Rattus norvegicus (Rat)
Target Names	TH

Image



Western Blot analysis of COLO205 cells using TH Polyclonal Antibody