





TRIM24 Antibody

Product Code	CSB-PA004291
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15164
Immunogen	Synthesized peptide derived from the C-terminal region of Human TIF1 α .
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TRIM24; RNF82; TIF1; TIF1A; Transcription intermediary factor 1-alpha; TIF1-alpha; E3 ubiquitin-protein ligase TRIM24; RING finger protein 82; Tripartite motif-containing protein 24
Immunogen Species	Homo sapiens (Human)
Target Names	TRIM24

Image



Western Blot analysis of COLO205 cells using TIF1 α Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using TIF1 α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.