




TNF Antibody

Product Code	CSB-PA004310
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01375
Immunogen	Synthesized peptide derived from the Internal region of Human TNF- α .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNF; TNFA; TNFSF2; Tumor necrosis factor; Cachectin; TNF-alpha; Tumor necrosis factor ligand superfamily member 2; TNF-a
Immunogen Species	Homo sapiens (Human)
Target Names	TNF

Image



Western Blot analysis of Jurkat cells using TNF- α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.