




TPH1 Antibody

Product Code	CSB-PA004319
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17752
Immunogen	Synthesized peptide derived from Human TPH1 around the non-phosphorylation site of S58.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TPH1; TPH; TPRH; TRPH; Tryptophan 5-hydroxylase 1; Tryptophan 5-monoxygenase 1
Immunogen Species	Homo sapiens (Human)
Target Names	TPH1

Image



Western Blot analysis of HepG2 cells using TPH1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.