




TRA2A Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA004320 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13595 |
| Immunogen | Synthesized peptide derived from the C-terminal region of Human Tra-2 α . |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | TRA2A; Transformer-2 protein homolog alpha; TRA-2 alpha; TRA2-alpha; Transformer-2 protein homolog A |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TRA2A |

Image



Western Blot analysis of HuvEc cells using Tra-2 α Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.