





TRADD Antibody

Product Code	CSB-PA004322
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15628
Immunogen	Synthesized peptide derived from the C-terminal region of Human TRADD.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TRADD; Tumor necrosis factor receptor type 1-associated DEATH domain protein; TNFR1-associated DEATH domain protein; TNFRSF1A-associated via death domain
Immunogen Species	Homo sapiens (Human)
Target Names	TRADD

Image



Western Blot analysis of COS7 cells using TRADD Polyclonal Antibody



Western Blot analysis of K562 cells using TRADD Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.