





NTRK2 Antibody

Product Code	CSB-PA004332
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16620
Immunogen	Synthesized peptide derived from Human Trk B around the non-phosphorylation site of Y706.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NTRK2; TRKB; BDNF/NT-3 growth factors receptor; GP145-TrkB; Trk-B; Neurotrophic tyrosine kinase receptor type 2; TrkB tyrosine kinase; Tropomyosin-related kinase B
Immunogen Species	Homo sapiens (Human)
Target Names	NTRK2

Image



Western Blot analysis of HepG2 cells using Trk B Polyclonal Antibody



Western Blot analysis of HepG2 cells using TrkB Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.