




TSSK4 Antibody

Product Code	CSB-PA004341
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6SA08
Immunogen	Synthesized peptide derived from the Internal region of Human TSSK 4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TSSK4; C14orf20; STK22E; TSSK5; Testis-specific serine/threonine-protein kinase 4; TSK-4; TSSK-4; Testis-specific kinase 4; Serine/threonine-protein kinase 22E
Immunogen Species	Homo sapiens (Human)
Target Names	TSSK4

Image



Western Blot analysis of HT29 cells using TSSK 4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.