




CACNA1H Antibody

Product Code	CSB-PA004342
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95180
Immunogen	Synthesized peptide derived from the Internal region of Human T-type Ca ⁺⁺ CP α1H.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CACNA1H; Voltage-dependent T-type calcium channel subunit alpha-1H; Low-voltage-activated calcium channel alpha1 3.2 subunit; Voltage-gated calcium channel subunit alpha Cav3.2
Immunogen Species	Homo sapiens (Human)
Target Names	CACNA1H

Image



Western Blot analysis of A549 cells using T-type Ca⁺⁺ CP α1H Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.