





TNFSF12 Antibody

Product Code	CSB-PA004351
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43508
Immunogen	Synthesized peptide derived from the Internal region of Human TWEAK.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFSF12; APO3L; DR3LG; Tumor necrosis factor ligand superfamily member 12; APO3 ligand; TNF-related weak inducer of apoptosis; TWEAK
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF12

Image



Western Blot analysis of 293 cells using TWEAK Polyclonal Antibody



Western Blot analysis of 293T AD293 HELA HepG2 KB cells using TWEAK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.