




UBA7 Antibody

Product Code	CSB-PA004359
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41226
Immunogen	Synthesized peptide derived from the C-terminal region of Human UBE1L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	UBA7; UBE1L; UBE2; Ubiquitin-like modifier-activating enzyme 7; Ubiquitin-activating enzyme 7; D8; Ubiquitin-activating enzyme E1 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	UBA7

Image



Western Blot analysis of HeLa cells using UBE1L Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.