





UBE3C Antibody

Product Code	CSB-PA004362
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15386
Immunogen	Synthesized peptide derived from the Internal region of Human UBE3C.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	UBE3C; KIAA0010; KIAA10; Ubiquitin-protein ligase E3C; HectH2
Immunogen Species	Homo sapiens (Human)
Target Names	UBE3C

Image



Western Blot analysis of HuvEc cells using UBE3C Polyclonal Antibody



Western Blot analysis of K562 cells using UBE3C Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.