




USP30 Antibody

Product Code	CSB-PA004395
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70CQ3
Immunogen	Synthesized peptide derived from the N-terminal region of Human USP30.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	USP30; Ubiquitin carboxyl-terminal hydrolase 30; Deubiquitinating enzyme 30; Ubiquitin thioesterase 30; Ubiquitin-specific-processing protease 30; Ub-specific protease 30
Immunogen Species	Homo sapiens (Human)
Target Names	USP30

Image



Western Blot analysis of HeLa cells using USP30 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.