




USP38 Antibody

Product Code	CSB-PA004396
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NB14
Immunogen	Synthesized peptide derived from the Internal region of Human USP38.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	USP38; KIAA1891; Ubiquitin carboxyl-terminal hydrolase 38; Deubiquitinating enzyme 38; HP43.8KD; Ubiquitin thioesterase 38; Ubiquitin-specific-processing protease 38
Immunogen Species	Homo sapiens (Human)
Target Names	USP38

Image



Western Blot analysis of 293 cells using USP38 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.