




USP50 Antibody

Product Code	CSB-PA004408
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70EL3
Immunogen	Synthesized peptide derived from the C-terminal region of Human USP50.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	USP50; Inactive ubiquitin carboxyl-terminal hydrolase 50; Inactive ubiquitin-specific peptidase 50
Immunogen Species	Homo sapiens (Human)
Target Names	USP50

Image



Western Blot analysis of A549 cells using USP50 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.