




UTP14A Antibody

Product Code	CSB-PA004409
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BVJ6
Immunogen	Synthesized peptide derived from the Internal region of Human UTP14A.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	UTP14A; SDCCAG16; U3 small nucleolar RNA-associated protein 14 homolog A; Antigen NY-CO-16; Serologically defined colon cancer antigen 16
Immunogen Species	Homo sapiens (Human)
Target Names	UTP14A

Image



Western Blot analysis of HeLa cells using UTP14A Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.