




VAV2 Antibody

Product Code	CSB-PA004442
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52735
Immunogen	Synthesized peptide derived from Human Vav2 around the non-phosphorylation site of Y142.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VAV2; Guanine nucleotide exchange factor VAV2; VAV-2
Immunogen Species	Homo sapiens (Human)
Target Names	VAV2

Image



Western Blot analysis of HeLa cells using Vav2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.