




VDR Antibody

Product Code	CSB-PA004443
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11473
Immunogen	Synthesized peptide derived from Human VDR around the non-phosphorylation site of S208.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VDR; NR111; Vitamin D3 receptor; VDR; 1; 25-dihydroxyvitamin D3 receptor; Nuclear receptor subfamily 1 group I member 1
Immunogen Species	Homo sapiens (Human)
Target Names	VDR

Image



Western Blot analysis of Jurkat cells using VDR Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.