




TNFSF15 Antibody

Product Code	CSB-PA004447
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95150
Immunogen	Synthesized peptide derived from the C-terminal region of Human VEGI.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFSF15; TL1; VEGI; Tumor necrosis factor ligand superfamily member 15; TNF ligand-related molecule 1; Vascular endothelial cell growth inhibitor
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF15

Image



Western Blot analysis of COLO205 cells using VEGI Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.