





VIM Antibody

Product Code	CSB-PA004449
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08670
Immunogen	Synthesized peptide derived from Human Vimentin around the non-phosphorylation site of S56.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VIM; Vimentin
Immunogen Species	Homo sapiens (Human)
Target Names	VIM

Image



Western Blot analysis of A549 cells using Vimentin Polyclonal Antibody



Western Blot analysis of 22RV1 HELA cells using Vimentin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.