





WNT1 Antibody

Product Code	CSB-PA004494
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04628
Immunogen	Synthesized peptide derived from the C-terminal region of Human Wnt-1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	WNT1; INT1; Proto-oncogene Wnt-1; Proto-oncogene Int-1 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	WNT1

Image



Western Blot analysis of HT29 cells using Wnt-1 Polyclonal Antibody



Western Blot analysis of NIH-3T3 HT29 cells
using Wnt-1 Polyclonal Antibody