





XIAP Antibody

Product Code	CSB-PA004519
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98170
Immunogen	Synthesized peptide derived from Human XIAP around the non-phosphorylation site of S87.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	XIAP; API3; BIRC4; IAP3; E3 ubiquitin-protein ligase XIAP; Baculoviral IAP repeat-containing protein 4; IAP-like protein; ILP; hILP; Inhibitor of apoptosis protein 3; IAP-3; hIAP-3; hIAP3; X-linked inhibitor of apoptosis protein; X-linked I
Immunogen Species	Homo sapiens (Human)
Target Names	XIAP

Image



Western Blot analysis of 293 cells using XIAP Polyclonal Antibody



Western Blot analysis of 4T1 293 cells using XIAP Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.