




ERCC5 Antibody

Product Code	CSB-PA004525
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28715
Immunogen	Synthesized peptide derived from the N-terminal region of Human XPG.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ERCC5; ERCM2; XPG; XPGC; DNA repair protein complementing XP-G cells; DNA excision repair protein ERCC-5; Xeroderma pigmentosum group G-complementing protein
Immunogen Species	Homo sapiens (Human)
Target Names	ERCC5

Image



Western Blot analysis of COLO205 cells using XPG Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.