




XRCC4 Antibody

Product Code	CSB-PA004533
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13426
Immunogen	Synthesized peptide derived from the C-terminal region of Human XRCC4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	XRCC4; DNA repair protein XRCC4; X-ray repair cross-complementing protein 4
Immunogen Species	Homo sapiens (Human)
Target Names	XRCC4

Image



Western Blot analysis of COLO205 cells using XRCC4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.