




ZBTB40 Antibody

Product Code	CSB-PA004562
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUA8
Immunogen	Synthesized peptide derived from the C-terminal region of Human ZBTB40.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZBTB40; KIAA0478; Zinc finger and BTB domain-containing protein 40
Immunogen Species	Homo sapiens (Human)
Target Names	ZBTB40

Image



Western Blot analysis of 293 cells using ZBTB40 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.