




ZC3H13 Antibody

Product Code	CSB-PA004565
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T200
Immunogen	Synthesized peptide derived from the C-terminal region of Human ZC3H13.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZC3H13; KIAA0853; Zinc finger CCCH domain-containing protein 13
Immunogen Species	Homo sapiens (Human)
Target Names	ZC3H13

Image



Western Blot analysis of A549 cells using ZC3H13 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.