





ZNF106 Antibody

Product Code	CSB-PA004567
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2Y7
Immunogen	Synthesized peptide derived from the C-terminal region of Human ZFP106.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZFP106; SH3BP3; ZNF474; Zinc finger protein 106 homolog; Zfp-106; Zinc finger protein 474
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF106

Image



Western Blot analysis of COS7 cells using ZFP106 Polyclonal Antibody



Western Blot analysis of COS-7 cells using ZFP106 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.