




ZRANB2 Antibody

Product Code	CSB-PA004584
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95218
Immunogen	Synthesized peptide derived from the Internal region of Human ZIS.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZRANB2; ZIS; ZNF265; Zinc finger Ran-binding domain-containing protein 2; Zinc finger protein 265; Zinc finger; splicing
Immunogen Species	Homo sapiens (Human)
Target Names	ZRANB2

Image



Western Blot analysis of HeLa cells using ZIS Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.