




ZNF134 Antibody

Product Code	CSB-PA004590
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52741
Immunogen	Synthesized peptide derived from the Internal region of Human ZNF134.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF134; Zinc finger protein 134
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF134

Image



Western Blot analysis of K562 cells using ZNF134 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.