




ZNF225 Antibody

Product Code	CSB-PA004591
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK10
Immunogen	Synthesized peptide derived from the C-terminal region of Human ZNF225.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF225; Zinc finger protein 225
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF225

Image



Western Blot analysis of HeLa cells using ZNF225 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.