




ZNF397 Antibody

Product Code	CSB-PA004594
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NF99
Immunogen	Synthesized peptide derived from the N-terminal region of Human ZNF397.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF397; ZNF47; ZSCAN15; Zinc finger protein 397; Zinc finger and SCAN domain-containing protein 15; Zinc finger protein 47
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF397

Image



Western Blot analysis of HuvEc cells using ZNF397 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.