




ZSCAN25 Antibody

Product Code	CSB-PA004604
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSZ9
Immunogen	Synthesized peptide derived from the C-terminal region of Human ZNF498.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF498; ZSCAN25; Zinc finger protein 498; Zinc finger and SCAN domain-containing protein 25
Immunogen Species	Homo sapiens (Human)
Target Names	ZSCAN25

Image



Western Blot analysis of COS7 cells using ZNF498 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.