




ZNF75D Antibody

Product Code	CSB-PA004610
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51815
Immunogen	Synthesized peptide derived from the C-terminal region of Human ZNF75.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF75D; ZNF75; ZNF82; Zinc finger protein 75D; Zinc finger protein 75; Zinc finger protein 82
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF75D

Image



Western Blot analysis of HuvEc cells using ZNF75 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.