




B3GALT1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA004621 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y5Z6 |
| Immunogen | Synthesized peptide derived from the Internal region of Human β -1,3-Gal-T1. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | B3GALT1; Beta-1; 3-galactosyltransferase 1; Beta-1,3-GalTase 1; Beta3Gal-T1; Beta3GalT1; UDP-galactose:beta-N-acetyl-glucosamine-beta-1,3-galactosyltransferase 1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | B3GALT1 |

Image



Western Blot analysis of HepG2 cells using β -1,3-Gal-T1 Polyclonal Antibody