




B3GALT4 Antibody

Product Code	CSB-PA004622
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96024
Immunogen	Synthesized peptide derived from the Internal region of Human β -1,3-Gal-T4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	B3GALT4; GALT4; Beta-1; 3-galactosyltransferase 4; Beta-1,3-GalTase 4; Beta3Gal-T4; Beta3GalT4; GalT4; b3Gal-T4; Gal-T2; Ganglioside galactosyltransferase; UDP-galactose:beta-N-acetyl-galactosamine-beta-1,3-galactosyltransferase
Immunogen Species	Homo sapiens (Human)
Target Names	B3GALT4

Image



Western Blot analysis of NIH-3T3 cells using β -1,3-Gal-T4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.