




B3GALTL Antibody

Product Code	CSB-PA004623
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6Y288
Immunogen	Synthesized peptide derived from the C-terminal region of Human β -1,3-Gal-TL.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	B3GALTL; B3GTL; Beta-1; 3-glycosyltransferase; Beta3Glc-T; Beta-3-glycosyltransferase-like
Immunogen Species	Homo sapiens (Human)
Target Names	B3GLCT

Image



Western Blot analysis of HepG2 cells using β -1,3-Gal-TL Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.