




ACADM Antibody

Product Code	CSB-PA004628
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11310
Immunogen	Synthesized peptide derived from the Internal region of Human MCAD.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACADM; Medium-chain specific acyl-CoA dehydrogenase, mitochondrial; MCAD
Immunogen Species	Homo sapiens (Human)
Target Names	ACADM

Image



Western Blot analysis of A549 cells using MCAD Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.