




CHD1L Antibody

Product Code	CSB-PA004631
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WJ1
Immunogen	Synthesized peptide derived from the Internal region of Human CHD1L.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CHD1L; ALC1; Chromodomain-helicase-DNA-binding protein 1-like; Amplified in liver cancer protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	CHD1L

Image



Western Blot analysis of NIH-3T3 Jurkat A549 cells using CHD1L Polyclonal Antibody