




FTO Antibody

Product Code	CSB-PA004640
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0B1
Immunogen	Synthesized peptide derived from the N-terminal region of Human FTO.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FTO; KIAA1752; Alpha-ketoglutarate-dependent dioxygenase FTO; Fat mass and obesity-associated protein
Immunogen Species	Homo sapiens (Human)
Target Names	FTO

Image



Western Blot analysis of K562 cells using FTO Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.