




ING4 Antibody

Product Code	CSB-PA004644
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNL4
Immunogen	Synthesized peptide derived from the Internal region of Human ING4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ING4; My036; Inhibitor of growth protein 4; p29ING4
Immunogen Species	Homo sapiens (Human)
Target Names	ING4

Image



Western Blot analysis of Jurkat K562 cells using ING4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.