




TRPV4 Antibody

Product Code	CSB-PA004651
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBA0
Immunogen	Synthesized peptide derived from the Internal region of Human TRPV4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TRPV4; VRL2; VROAC; Transient receptor potential cation channel subfamily V member 4; TrpV4; Osm-9-like TRP channel 4; OTRPC4; Transient receptor potential protein 12; TRP12;
Immunogen Species	Homo sapiens (Human)
Target Names	TRPV4

Image



Western Blot analysis of K562 cells using TRPV4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.