




GPM6A Antibody

Product Code	CSB-PA004657
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51674
Immunogen	Synthesized peptide derived from the Internal region of Human M6A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPM6A; M6A; Neuronal membrane glycoprotein M6-a; M6a
Immunogen Species	Homo sapiens (Human)
Target Names	GPM6A

Image



Western Blot analysis of rat stomach cells using M6A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.