




STAT2 Antibody

Product Code	CSB-PA004664
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52630
Immunogen	Synthesized peptide derived from Human Stat2 around the non-phosphorylation site of Y631.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STAT2; Signal transducer and activator of transcription 2; p113
Immunogen Species	Homo sapiens (Human)
Target Names	STAT2

Image



Western Blot analysis of K562 cells using Stat2 Polyclonal Antibody