




ACTA2 Antibody

Product Code	CSB-PA004667
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62736
Immunogen	Synthesized peptide derived from the N-terminal region of Human α -SMA.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACTA2; ACTSA; ACTVS; GIG46; Actin; aortic smooth muscle; Alpha-actin-2; Cell growth-inhibiting gene 46 protein
Immunogen Species	Homo sapiens (Human)
Target Names	ACTA2

Image



Western Blot analysis of Jurkat K562 HuvEc cells using α -SMA Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.