




APOL6 Antibody

Product Code	CSB-PA004678
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWW8
Immunogen	Synthesized peptide derived from the C-terminal region of Human ApoL6.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APOL6; Apolipoprotein L6; Apolipoprotein L-VI; ApoL-VI
Immunogen Species	Homo sapiens (Human)
Target Names	APOL6

Image



Western Blot analysis of K562 cells using ApoL6 Polyclonal Antibody