




FAIM2 Antibody

Product Code	CSB-PA004691
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWQ8
Immunogen	Synthesized peptide derived from the N-terminal region of Human FAIM2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FAIM2; KIAA0950; LFG; LFG2; NMP35; TMBIM2; Protein lifeguard 2; Fas apoptotic inhibitory molecule 2; Neural membrane protein 35; Transmembrane BAX inhibitor motif-containing protein 2
Immunogen Species	Homo sapiens (Human)
Target Names	FAIM2

Image



Western Blot analysis of Jurkat cells using FAIM2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.